

PVMS OF MEDICAL EQUIPMENT

Clinical Specialty

Radiological Equipment

Generic Name

X-RAY FILM PROCESSOR

Clinical Purpose

It helps in automatic processing of x-ray films providing better quality images and consuming less time

TECHNICAL SPECIFICATIONS

High throughput free Standing type Continuous roller transport Sheet film,
10x10cm-35x43cm sizes (14"x17")

High through put more than 110 sheets (35 x 43cm)

Forced hot air circulation method

Piping installed in the DEV and FIX racks.

Processing solution temperature; Controlled by the temperature control tank, with the thermistor monitoring, with the heater heating and with the wash water cooling.

Drying temperature; Controlled automatically according to a fixed temperature setting.